

Sung Kwang Meditech

For surgical patient & Rehabilitation Radiology Products

Operating Room


Snow pac (Bean Bag) Several device Position pad Gel product Instrument Case Silicon product



Rehabilitation Gel wheelchair cushion pad Gel wheelchair back cushion pad

Gel embo cushion mattress





Radiation Protection For: Cardiac Cath Lab Interventional Radiology Pain Management









T) 82-31-855-7548 F) 82-31-855-7549 E) skmeditech@skmeditech.com A) 1014-20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea




SNOW PAC (Bean Bag)

FOR SURGERY POSITIONING PATIENTS

Snow Pac are soft pads (Available in various sizes and shapes) that are filed with thousands of tiny plastic beads. Snow Pac provides more stable and comfortable position to patients and saves valuable time for the entire O.R team on every case.

It aids to become comfortable position. It aids better surgery position presentation.

PRINCIPLE OF SYSTEM





EFFECT

Stable

Because Snow-Pacs are molded well, firmly and securely to shape of the patient's body than sandbags, towel rolls and other devices that are merely propped against the patient.

Comfortable

The apparent wrinkles are inverted-drawn inward by suction and will not irritate the skin.

Better Presentation

With Snow-Pacs the patiant can be positioned exactly as the surgeon desires for excellent presentation of the operative site. Snow-Pacs conform to any shape, so the patient's position does not have to be compromised.

More Versatile

With Snow Pacs patient can often be placed in positions that could not be easily achieved using other positioning devices.

No Pressure Points (With Gels)

A Snow-Pac supports all parts of the patient's body uniformly when used with gel like a floating mattress.

Saves Time

Patients can be positioned much more quickly and easily using Snow Pacs, and if necessary, be repositioned during a case.

Before

Fix or put the Snow Pac to the patient's body according to the position you want to form (e.g., Supine Position) and connect the suction line to the entrance of the Snow Pac.

Depending on the intensity of the Snow Pac you want to make, use a suction line to inhale the air in the Snow Pac properly. (If you inhale for about 10 to 15 seconds, vaccum is created inside the Snow Pac and it becomes hardens)

After

The inside of the Snow Pac is vacuumed and fixed to the desired shape.

When the suction line is disconnected, the Snow-Pac retains its firm, molded shape and holds the patient securely in position.

After surgery, if user open the entrance of the Snow Pac and press it slightly, air re-enters into the Snow Pac and it becomes soft again and re-useable.

OTHER FEATURES

Any Vacuum Source

Use wall suction, aspirator pump or any other convenient O.R. vacuum source.

Reusable

As the exterior material is soft and sturdy, Snow Pac can be used for many years as long as they are thoroughly controlled for contamination.

Radiolucent

Excellent radio transparency with only minimal and clearly identifiable artifacts.

Clean, Disinfect, Sterilize

Can be cleaned with regular soaps and liquid disinfectants or sterilized by gas.

Saves Money

Snow Pacs save the cost of sandbags, towel rolls and other positioning devices, and they save time in every case.

Repair Service

Factory inspection and repair service is available at nominal cost after warranty period. Full One Year Warranty.

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

INDICATION

KIDNEY POSITION, LATERAL POSITION





A-1 type A-1-1 type



A-1-2 type A-1-3 type

Easily up in nephrectomy (A type), and also for handicapped person with cerebral palsy.(A-1 type).



A type	73(L) * 100(W)cm
A-1 type	73(L) * 100(W)cm
A-1-1 type	90(L) * 100(W)cm
A-1-2 type	73(L) * 100(W)cm(U)
A-1-3 type	90(L) * 100(W)cm(U)

Provides hyperextension, prevents postoperative discomfort.

B-1 type 118(L) * 88(W)cm	

THYROID POSITION



B type B-1 type

HEAD STABILIZATION POSITION



SUPINE POSITION



Stabilizes head in eye, facial or dental surgery.


C type 30(L) * 40(W)cm

Supports and eliminates back pains during long procedures.

Prevents foot drop splint in Emergency.


D type 30(L) * 55(W)cm





D type

POSITION PADS



Tel: 82-31-855-7548 Fax: 82-31-855-7549

Autoclave Instrument Cases

Instrument Sterilization Cases are designed to protect fine surgical instruments and various scopes during sterilization and storage. Our cases are made of elastic special materials for flexible handling and **it is transparent so that the inside can be seen during autoclaving.** There are fine circular holes in the upper(Lid) and lower(Base) of the case to **help circulate steam during sterilization.** Silicon Embo Mat or Bar is placed on the bottom and prevents instruments from moving freely.


- This sterilizing instrument tray is designed to withstand repeated exposure to 134 degrees (Steam / EO / Plasma steralizable).
- There are a variety of designs and sizes suitable for device applications such as **micro instrument**, Laparoscope, Telescope, and Endoscope.
- The size table below is based on the outer diameter, and the inner diameter by -1 to 2cm.

Micro Instrument	Case(26(L) x 17(W))	Scope Case(Telesc	ope, Laparoscope)
SK-19(2F)	26(L) x 17(W) x 4(H) cm	SK-20 (Long Scope Case)	60(L) x 10(W) x 5(H) cm
SK-19-a(1F)	26(L) x 17(W) x 2.5(H) cm	SK-20-1 (Laparoscope Case)	48(L) x 9(W) x 5(H) cm
		SK-20-2(Telescope Case)	28(L) x 7(W) x 5(H) cm
	К-19(2F)	SK-20-1	
Micro Instrument	Case(1F)	Micro Instrument C	Case(1F, 2F)
SK-19-S	16(L) × 6(W) × 3(H) cm	SK-19-5(2F)	39(L) x 27(W) x 4(H) cm
SK-19-1	25(L) x 5(W) x 2.5(H) cm	SK-19-5-a(1F)	38(L) x 26(W) x 2.5(H) cm
SK-19-3	30(L) x 10(W) x 2.5(H) cm	SK-19-5-b(1F)	39(L) x 27(W) x 12(H) cm
	SK-19-S	SK-19-5-a	SK-19-5-b
Scope Case(Long,	Big size)	Ophthalmic Instrumen	t Case(Silicon color inside: blue)
SK-20-3(Endoscope)	37(L) x 15(W) x 7(H) cm	SK-19-a(1F)	26(L) x 17(W) x 2.5(H) cm
SK-20-4(Laparoscopic)	$60(L) \times 20(W) \times 6(H) \text{ cm}$	SK-19-S	16(L) x 6(W) x 3(H) cm
SK-20-5(Long scope)	72(L) × 15(W) × 10(H) cm	SK-19-5-a(1F)	38(L) × 26(W) × 2.5(H) cm
SK-21-1 (Laparoscopic)	64(L) × 25(W) × 15.5(H) cm		
SK-21-1-A (Laparoscopic)	64(L) x 25(W) x 12.5(H) cm	Micro Instrument Case	K-19-a
Sung Kwang Meditech Co.,	Ltd www.skmeditec	h.com Contact to	E-mail: skmeditech@skmeditech.com
			-855-7548 Fax: 82-31-855-7549
			5

CASES

INSTRUMENT STERILIZATION TRAYS USE CASES

The instrument sterilization case photos below are examples of the use of cases that were actually applied in storage or sterilization processes by various departments such as surgery and internal medicine.




SK-19-5(2F) 39(L) x 27(W) x 4(H) cm



SK-20-3 37(L) x 15(W) x 7(H) cm



SK-19-a(1F) 26(L) x 17(W) x 2.5(H) cm



SK-19-5-b(1F) 39(L) x 27(W) x 12(H) cm

Sung Kwang Meditech Co.,Ltd



SK-21-1 64(L) x 25(W) x 15.5(H) cm



SK-19-3 30(L) x 10(W) x 2.5(H) cm



Strengths of Sungkwang Meditech trays:

- Cases are transparent. (Users can see instruments during autoclaving)
- Cases are light but strong. (We do not use stainless, steel materials. But, use special strong materials)
- There are fine circular holes in the upper(Lid) and lower(Base). (It helps circulate steam during sterilization)
- Silicon Embo Mat or Bar is available. (It is placed on the bottom and prevents instruments from moving freely).
- It is possible to assemble freely based on the shape of the instrument.



SK-21-1 64(L) x 25(W) x 15.5(H) cm



SK-20-3 37(L) x 15(W) x 7(H) cm



SK-20-5 72(L) x 15(W) x 10(H) cm



SK-20-4 60(L) x 20(W) x 6(H) cm



SK-20-4 60(L) x 20(W) x 6(H) cm



SK-20-3 37(L) x 15(W) x 7(H) cm



SK-20-3 37(L) x 15(W) x 7(H) cm



SK-20-3 37(L) x 15(W) x 7(H) cm



SK-20-4 60(L) x 20(W) x 6(H) cm

CASES

NEEDLE COUNTER & SHARPS COLLECTOR

Inner : magnet For surgery in OR It aids to count needle thoroughly.



NEEDLE HUNTER It finds needles by rolling on the floor



		(External)
SK-18 (S)	14(L) * 4(W) * 1.5(H)cm	
SK-18-1(L)	17(L) * 8(W) * 2(H)cm	
SK-13-2	100(L) * 30(W)	

MICRO INSTRUMENT ONE FLOOR CASE

SK-18





(External) SK-19-a(1F) 26(L) * 17(W) * 2.5(H)cm

MICRO INSTRUMENT TWO FLOOR CASE



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposureto 134°C steam temperature.

Steam / EO / Plasma Sterilizable

SK-19(2F) 26(L) * 17(W) * 4(H)cm (External)

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

MICRO INSTRUMENT ONE FLOOR CASE



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.



Steam / EO / Plasma Sterilizable



(External) SK-19-S 16(L) * 6(W) * 3(H)cm SK-19-1 25(L) * 5(W) * 2.5(H)cm SK-19-3 30(L) * 10(W) * 2.5(H)cm

Sung Kwang Meditech Co.,Ltd

Contact to E-mail: skmeditech@skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549

MICRO INSTRUMENT ONE FLOOR CASE



MICRO INSTRUMENT TWO FLOOR CASE

For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument

Cases are designed to withstand repeated exposure to 134°C steam

Steam / EO / Plasma Sterilizable

(External)

39(L) * 27(W) * 4(H)cm


to move freely. temperature. SK-19-5 (2F) SK-19-5(2F) Contact to E-mail: skmeditech@skmeditech.com

```
Sung Kwang Meditech Co.,Ltd
```

www.skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549





Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549

For Telescope.

ENDOSCOPE INSTRUMENT CASE



SK-20-3





For Endoscope Instrument.

Silicon bar is placed to prevent laparoscopic instrument to move freely.

Cases are designed to with stand repeated exposure to 134°C steam temperature.

SK-20-3 is used for endoscope camera case.

Steam / EO / Plasma Sterilizable



SK-20-3 37(L) * 15(W) * 7(H) cm

Silicon bar is placed to prevent laparoscopic instrument to move

Cases are designed to withstand repeated exposure to 134°C steam

Steam / EO / Plasma Sterilizable

It contains 2pcs or more of laparoscopic instrument.

freely.

temperature.

(External)

(External)

LAPAROSCOPIC INSTRUMENT CASE



SK-20-4


SK-20-4 60(L) * 20(W) * 6(H)cm

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

LONG SCOPE INSTRUMENT CASE





LAPAROSCOPIC INSTRUMENT CASE

Sterilization cases used during laparoscopic long scope and robotic instrument surgery.

Silicon bar or mat is placed to prevent movement of instruments used in obstetrics, urology, otorhinolaryngology, surgery and robotic surgery.

Two or more laparoscopic instruments can be inserted.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)


SK-20-5 72(L

Tel: 82-31-855-7548 Fax: 82-31-855-7549

72(L) * 15(W) * 10(H)cm



OPHTHALMIC INSTRUMENT CASE



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

		(External)
SK-19-S	16(L) * 6(W) * 3(H)cm	
SK-19-a (1F)	26(L) * 17(W) * 2.5(H)cm	
SK-19-5-a (1F)	38(L) * 26(W) * 2.5(H)cm	

BIOPSY CASE

-			
R2	Delphain	Pyramidal	L2
R3			L3
R. Jungular	R. Thyroid	L. Thyroid	L. Jungular
R4			L4
R5	R.CN	L.CN	L5
R. Para. T			L. Para. T

SK-21-a-1 (for Thyroid)

Rt Le	vel 1	Rt Level 1					
Rt Level 2		Rt Level 2					
Rt Level 3		Rt Level 2					
Rt deep margin	Rt deep margin Rt rotter's		Lt rotter's				
Rt. axillary	Rt. axillary Rt. sentinal		Lt. sentinal				
SK-21-2-2 (for broast)							

SK-21-a-2 (for breast)

	22	Delphian	Pyramidal	L2
	83 2721 84	A.Thyroid	L. Thyroid	L3
	RS	R.CN	L.CN	1.5
-	- an			Pare

	Rtie	vel 1	Lt level 1	
• •	Rtie	vel 2	LLL JOY	/012
	Rt ie	yel 3	LUNE	13
1	na deep minraîn	Heretter s	deep mirin	conters.
- 0	REatling	lenina a-	Wallian	Internet.

Material : stainless steel with marked
L/N number,
Use : for specimen storage

·

order made available

(External)


 SK-21-a-1
 for Thyroid 16.5(W) * 19(L) * 3.2(H)

 SK-21-a-2
 for breast 14.5(W) * 15(L) * 3.2(H)

 SK-21-a-3
 for stomach 31.5(W) * 8.5(L) * 3.2(H)

 SK-21-a-4
 for biliary and Pancreas 30(W) * 8.5(L) * 3.2(H)

	1	2	3	4	5	6	7	4sa	1 a b	44
	8	9	10	11	11p	11d	12	14v	450	40
SK-21-a-3 (for stomach)										



SK-21-a-4 (for biliary and Pancreas)



Contact to E-mail: skmeditech@skmeditech.com

www.skmeditech.com

GEL PRODUCTS

Capillary closure is at the root of ulcer formation.

The blockage or damage of capillaries prevents the natural maintenance of healthy skin.

In operation, the use of proper positioning devices prevent excessive pressure build-up and capillary closure.

Evenly distribute weight and provide constant support beneath patients.

Gel products are cost effective, latex free, fast & easy cleaning.




Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

Sung Kwang Meditec	h Co.,Ltd wv	ww.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.co Tel: 82-31-855-7548 Fax: 82-31-855-75
PRONE GEL HEAD P	SK-22-2	SK-22-2-1	It provides clear air passage way when the patient is positioned for surgery by prone position. SK-22-2 27(L) * 24(W) * 15/12(H)cm SK-22-2-1 27(L) * 21(W) * 8(H)cm
	JV-77_6	SK-22-5	SK-22-e(Adult) 18(L) * 18(W) * 8(H)cm SK-22-5(Baby) 16(L) * 12(W) * 1(H)cm
GEL HEAD PAD	SK-22-e		It provides optimum stability of patient s head while in the supine position. Can also be used for position and protect the knee or calf.
	SK-22-c	SK-22-d	SK-22-c(Adult)23(L) * 23(W) * 5.7(H)cmSK-22-d(Child)19(L) * 21.6(W) * 3.8(H)cm
			It helps prevent risk of alopecia especially during long surgery. Cradling of patient s head provides stability.
GEL OPHTHALMIC CRADLE HEADREST			It was designed for ENT, plastic surgery ,all ophthalmic cases.
	SK-22-a	SK-22-b	SK-22-a(Adult)31(L) * 27(W) * 10(H)cmSK-22-b(Child)31(L) * 25(W) * 7.3(H)cm
			Also provides patient with comfort while under local anesthesia, minimizing patient movement. It was also designed for surgeon's functional arm support.
GEL OPHTHALMIC HEADREST			It was designed for stable head placement for all ophthalmic and ora or facial procedures.
	SK-22-1	SK-22-1-2 SK-22-1-1	SK-22-1 (Adult) 21(D) * 4.5(H)cm SK-22-1-1 (Child) 15(D) * 3.5(H)cm SK-22-1-2 (Baby) 9.2(D) * 2.4(H)cm
HEAD PAD		60	clear air passage way when the patient is positioned for surgery by lateral or prone position.
C-DOUGHNUT GEL			C-doughnut gel head pad provides
	SK-22	SK-22-4 SK-22-3	SK-22-3(Child) 15(D) * 3.5(H)cm SK-22-4(Baby) 9.2(D) * 2.4(H)cm
			patient is positioned for surgery by supine position. SK-22(Adult) 20(D) * 4.5(H)cm
HEAD PAD			excellent protection when the

Sung Kwang Meditech $C \in$



Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

GEL KNEE PAD			surgery in	pad is used for OBGY lithotomy position. s against nerve injury.
	SK-25 (S) SK-25-1 (L)	all and the second s	SK-25 SK-25-1	40(L) * 24(W) * 1(H)cm 50(L) * 52(W) * 1(H)cm
GEL LATERAL PAD GEL				I pad is designed to
		GEL THYROID PAD	provides e	date both the shoulder and exceptional pressure relief e protection.
		T	thyroidec with gent	d pad is designed for tomy. It provides patient le pressure while long ith hyperextension position
	SK-26	SK-26-2	The tissue comfortal	e friendly feels make it very ole.
			SK-26 SK-26-2	70(L) * 50(W) * 12(H) / 0.8cm 45(L) * 30(W) * 7(H)cm
GEL PATIENT POSITIONER	SK-2	27-a SK-27-b	lt is good versatile p	for any procedure. positioner.
- Constant		SK-27-c	SK-27-a(Sm SK-27-b(Lar SK-27-c(Ext	ge) 35(L) * 12(W) * 3(H)cm
GEL PRONE CHEST PAD		SK-27-e	conform t	e designed for uniformly to a patient s shoulder to a, with rounded edges for
				ide optimum alleviating on the thorax area.
			SK-27-d(Rig SK-27-e(Lef	
GEL SHOULDER PAD (CHEST)			maximum	ular design provides pressure protection while in supine position.
				4 can be used under the foot against the board.
SK-27				ient is in lateral position, be used for the shoulder n.
		1.		ectomy , SK-27-4 is used fo erextension.
SK-27-1 (hard)	I		while pati	ent is in prone position SK- 7-2 can be used under the
	A. 1981		SK-27 SK-27-1	47(L) * 10(W) * 5(H)cm 47(L) * 15(W) * 8(H)cm
SK-27-3	SK-27-2 (hard)	SK-27-4	SK-27-2 SK-27-3	47(L) * 15(W) * 12(H)cm 25(L) * 15(W) * 8(H)cm
			SK-27-4	52(L) * 10(W) * 8(H)cm

Sung Kwang Meditech Co.,Ltd


www.skmeditech.com

Sung Kwang Meditech C \in



GEL FOREARM PROTECTOR		It is used for nerve protection of arm and elbow in supine position surgery.
		SK-17-2 Gel forearm protector (L)
		46(L) * 21(W) * 0.8(H)cm
SK-17-2 (L)	As a Asia	SK-17-2-1 Gel forearm protector (S)
SK-17-2-1 (S)		38(L) * 11.5(W) * 0.8(H)cm
GEL BODY BELT		It helps safety, stability and tissue protection during surgery.
SK-17-3		SK-17-3 50(L) * 8.5(W) * 0.8(H)cm
GEL SHOULDER PROTECTOR		It is designed for patient in trendelenburg down position to protect shoulder areas.
	SK-17-3-1	SK-17-3-1 16(L) * 15(W) * 2(H)cm
GEL ALLEN PAD SK-17-4-1		It is used lithotomy position for nerve protection.
		SK-17-4-1 68(L) * 38(W) * 1(H)cm
GEL ANKLE HOLDER (2 PIECE SET)		It is used in lithotomy position for nerve protection.
	SK-17-5	SK-17-5 40(L) * 3(W) * 1(H)cm
GEL SNOW PAC A-1 type COVER		Designed for A-1 type snow pac.
		It provides comfort.
	SK-17-6	SK-17-6 100(L) * 73(W) * 1(H)cm
GEL SNOW PAC B type COVER		Designed for B type snow pac.
		SK-17-6-1 88(L) * 72(W) * 1(H)cm
	SK-17-6-1	
GEL SNOW PAC C type COVER		Designed for°C type snow pac.
	SK-17-6-2	SK-17-6-2 30(L) * 40(W) * 1(H)cm
	JN-1/-0-Z	3n 1/ U 2 3U(L) ^ 4U(W) * 1(D)(III
Sung Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.com
20		Tel: 82-31-855-7548 Fax: 82-31-855-7549

SILICON



Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com

6.0(I.D) * 11.0(O.D)mm

21

Tel: 82-31-855-7548 Fax: 82-31-855-7549

SK-31-4-4

Several Devices

PATIENT TRANSPORT BOARD



The fully equipped patent transport board is ideal for transport when patients need to be transferred from MRI, operating table, bed, etc.

SIZE:	
SK-MTB80	80 x 35cm
SK-MTB140	140 x 45cm
SK-MTB183	183 x 55cm

RADIATION SHIELDING X-RAY CURTAINS FOR SCANNER



Radiation Shielding Curtains for security, food scanners. The protection range includes lead equivalent values of 0.125, 0.25, 0.375, 0.5mm depending on the level of protection required. Customer made.

ARM PROTECTOR





SK-CU-CP(S)

SK-13

SK-14

SK-14-1

SK-14-6

SK-14-7

Multi-Purpose position pillow for hospital operating room.

It is easy to use for short-term surgery like Hemorrhoid surgery

Inner : soft cotton cushion

40.5(L) * 17(H) * 15.5(W)cm Color : White

SK-CU-CP(L)

PROTECTIVE BODY BAND

PROTECTIVE KNEE BAND



Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

200(L) * 11(W)cm

18(L) * 18.5(W) cm

35(L) * 12(W) cm

70(L) * 7(W)cm

 SK-CU-CP-L
 45(L) * 70(W) cm


 SK-CU-CP-S
 40(L) * 60(W) cm

Custom made available

GEL WHEELCHAIR CUSHION PAD WITH COVER

It was made with a high-quality polyurethan gel and PU foam. The inner gel material is highly resistant to bacteria and bodily fluids, and can be easily cleaned with warm soapy water or a disinfectant. Also, it helps dissipate heat buildup.

Gel Wheelchair cushion pad offers excellent pressure reduction & distribuion, great comfort, and positioning stability. It can be applied for individuals with a history of skin breakdown or those who suffer from pressure sores. Feel free to enjoy our products below.



GEL WHEELCHAIR BACK CUSHIONS WITH COVER

Polyurethane gel and PU foam helps comfort and durability for long-term.

The embossed gel feature aids comfort for long-lasting back and lumbar support in long-term seating environments.



SNOW PAC



It is used for lugeric patient and unpositioned handicapped person. Make position easy and fast. Then, suction equipment is used.

RH B type	83(L) * 108(W)cm
HB-102a	33cm

GEL HORSE S EQUESTRIAN PAD

RECTANGULAR ROUND




SK-EQ-a Rectangular shape with cover
36(L) * 50(W) * 1(H) cm
SK-EQ-a cover Rectangular shape cover only
36(L) * 50(W) * 1(H) cm
SK-EQ-b Round shape with cover
40(L) * 55(W) * 1(H) cm
SK-EQ-b cover Round shape cover only
40(L) * 55(W) * 1(H) cm

SILVER PRODUTS

GEL EMBO CUSHION PILLOW / MATTRESS WITH COVER (Memory foam + embo gel)

Characteristics:

C-spine supporting effect is developed by medical basis. As humanengineering structure, it softly supports head and provides comfortable.



Outer : 100% cotton velvet with silvercholoid coating Inner : cotton 75% polyester 25%

nner Collon 75% polyester 25%

Materials : slowly recovery polyurethane foam 100% + embo gel pad

Color : silver, ivory

Kind : therapeutic, antibacterial, deodorizational.

RH-1	50(L) * 33(W)cm
RH-2 (twin)	RH-2-1 (full)
RH-2-3 (queen)	RH-2-4 (king)
Cover Only	
RH-2-cover	RH-2-1-cover
RH-2-3-cover	RH-2-4-cover

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

Lead Aprons

Several Types of Aprons For Radiation Protection



SK-15 (Buckle type)



SK-15-1 (Coat type)





SK-15-4 (X-band type)



SK-15-5 (Buckle X-band type)

Sung Kwang Meditech Co.,Ltd



Here are all the radiation protection products you need for protection from radiation exposures.

We manufacture a comprehensive range of high-quality, lightweight, and cost-effective products for medical customers around the

world.

Inner : Natural rubber-based Lead flexible wearing comfort and best.

Outer : nylon fabric, etc satisfies the workmanship hygienic requirements and the demands of medical use.

Pb Equivalent :

0.25 / 0.35 / 0.50 mmPb Available

±2~4cm, unit∶cm MALE L:Length,W:width Μ XL S L W W W W L L L 103 Buckle type 93 58 98 60 63 108 66 Coat type (Front/back) 55 105/95 96/86 53 99/89 57 106/96 60 X-band type 94 57 100 60 104 55 96 62 51/43 Vest/Skirt type 59/50 49/42 60/52 61/53 53/46 63/54 56/49 82 47 85 50 88 53 91 55 panorama type

FEMALE

	9	5	N	Л	L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.

Available to SK-15, SK-15-1, SK-15-4, SK-15-5

www.skmeditech.com

SK-15-6 (X-buckle type)

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

Lead Aprons

Several Types of Aprons For Radiation Protection



Tel: 82-31-855-7548 Fax: 82-31-855-7549 27

Non Lead Aprons

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15 (Buckle type)





SKN-15-1(Coat type)





SKN-15-4 (X-band type)



SKN-15-5 (Buckle X-band type)

Sung Kwang Meditech Co.,Ltd



Here are all the radiation protection Lead-Free products you need for protection from radiation exposures.

We manufacture a comprehensive range of high-quality, lightweight, and cost-effective products for medical customers around the world.

- Inner : Natural rubber-based Lead flexible wearing comfort and best.
- Outer : Water Resistant, Abrasion Resistant, Cool wear

Pb Equivalent :

0.25 / 0.35 / 0.50 mmPb Available

	world.			MA	ALE .			cm, unit∶cm th, W∶width
	9	5	N	Λ	L	_	X	L
	L	W	L	W	L	W	L	W
Buckle type	93	58	98	60	103	63	108	66
Coat type (Front/back)	96/86	53	99/89	55	105/95	57	106/96	60
X-band type	94	55	96	57	100	60	104	62
Vest/Skirt type	59/50	49/42	60/52	51/43	61/53	53/46	63/54	56/49
panorama type	82	47	85	50	88	53	91	55

FEMALE

	9	5	N	Л	L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.

Available to SKN-15, SKN-15-1, SKN-15-4, SKN-15-5

www.skmeditech.com

SKN-15-6 (X-buckle type)

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

Non Lead Aprons

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15-3 (Half type)

SKN-16



SKN-15-2 (Vest & Skirt type)

SKN-15-3-1 (Panorama type)

SKN-16-a

SKN-16-b



SKN-15-3-2

Various types of aprons available.

Thyroid shield is an essential to protect the thyroid gland from radiation exposure. (SKN-16, SKN-16-a, SKN-16-b)

Find a wide variety of radiation protection products and colors.

Pb Equivalent: 0.25 / 0.35 / 0.50 mmPb Available

Color type





Dia Chocolate



Dia Orange



Pearl Purple

Pearl Navv



Pearl Grev













Melange Navy





Black Dot

Honeycomb Grey

Our Lead-Free cap and gonad shield are lightweight and comfortable. An adjustable velcro in the back is attached for a customizable fit. -Two types for cap : SKN-16-4 / 16-4-a -Two types for gonad : SKN-16-5 / 16-5-a -Order Types : Male, female

-Protection : 0.25mm front

-Pb Equivalent: 0.25 / 0.35 / 0.50 mmPb Available

Custom-made available.

Lead Free cap SKN-16-4 (Basic Model) Lead Free cap SKN-16-4-a (Cool wear mositure-wicking fabric) Lead Free gonad shield SKN-16-5 Lead Free gonad shield SKN-16-5-a

Sung Kwang Meditech Co.,Ltd

SKN-16-5

SKN-16-4

www.skmeditech.com

SKN-16-5-a

SKN-16-4-a



APRON HOLDER

APRON WALL CONSOLE



Apron holder will accommodate up to 2-10 aprons.

Unit is ideal if wall space is limited.

It is designed to prevent declining on side in moving.

Solid stainless steel construction prevents rust in long-term use.

High durability under heavy use.

Can be assembled after delivering.

SK-16-6	up to 10
SK-16-6-a	up to 5
SK-16-6-1	up to 2
SK-16-6-2	up to 5

FOR PREGNANT WOMAN BARRIER TABLE LEAD / NON-LEAD CURTAIN







SK-16-7-4 (양면) SK-16-7-4 -a (단면) Lead curtain provides separated radiation shielded space by simply drawing.

- protection : 0.5 mmpb
- Order made
- Protection for pregnant woman barrier : 0.5 mmpb
- Order made

SK-16-6	71(L) * 63(W) * 140(H) cm
SK-16-6-a	50(L) * 63(W) * 140(H) cm
SK-16-6-1	20(L) * 50(W) * 25(H) cm
SK-16-6-2	25(L) * 50(W) * 25(H) cm
SK-16-7-1	57(W) * 60(L)/83(H) to 146(H) cm
SK-16-7-3	200(W) * 150(H) cm
SK-16-7-4	100(L) * 65(W) * 70(H) cm
SK-16-7-4-a	100(L) * 65(W) * 70(H) cm

LEAD BARRIER



Lead-glass viewing panel provides exceptional visibility.

- Protection : 1.76 mmpb
- Order made

SK-16-7 lead glass barrier	180(L) * 90(W) cm
Window Size:	80(L) * 90(W) cm
SK-16-7-a lead barrier	180(L) * 90(W) cm
SK-16-7-b	180(L) * 90(W) cm
Window Size:	20(L) * 30(W) cm
SK-16-7-c	180(L) * 90(W) cm
Window Size:	180(L) * 90(W) cm

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Radiation Protection Eyewear

Since **human's eyes are more sensitive** to radiation than other parts of human body, **using eyewears or masks is strongly recommended** for users related to the following vocations: Fluoroscopic Procedures, Cardiac Cath, EP Labs, Urology, Interventional Radiology, Orthopedic Surgery.

- Lead Masks for Full Faces: 0.1 mm Pb
- Lead Eyewears (With Glasses / Without Glasses): 0.75 mm Pb

Lead Masks

SK-400PSM (PANORAMIC FACE SHIELD)



The 400PSM Panorama Shield consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from a dual adjustable headpiece.

The acrylic shield can be positioned at any angle to the face so as to minimize scatter radiation from passing under the shield by side adjustable dials. The headgear is adjustable in two positions.

An adjustable headband that saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back. The circumference portion of the headband has an adjustable dial in the rear to increase or decrease the hat size for a comfortable, snug fit.

Wipe clean with a soft cloth. Abrasive materials may scratch front lens. Cannot be sterilized with high temperatures or high pressure.

- ◆ Weight : 350 grams
- ◆ Protection : .10mm Lead Acrylic
- ◆ Material : Acrylic/Nylon

SK-450FPM (FULL FACE SHIELD)



The 450FPM Panorama Mask consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from an adjustable headpiece.

The stationary acrylic shield extends from a headpiece brow connector for a length of 8 inches from top to bottom in front of the entire face. A support strap saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back.

The circumference portion of the headband has an adjustable strap in the rear to increase or decrease the hat size for a comfortable, snug fit. W ipe clean with a soft cloth.

Abrasive materials may scratch front lens. Cannot be sterilized with high temperatures or high pressure.

- ◆ Weight: 450 grams
- ◆ Protection: .10mm Lead Acrylic
- ◆ Material: Acrylic/Nylon

SK-450FPM

SK-RE-FGUARD (Beta/Alpha)



The Fitguard fitover is a versatile, high guality nylon frame. Featuring Flexfit technology, the Fitguard lead glasses' temples can be benched (adjusted) easily for a secure, comfortable fit. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

The Beta model is ideal for fitting over standard prescription eyewear with square or rounded lenses.

The Alpha model is ideal for fitting over wide prescription evewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the Fitguard Alpha is our most popular fitover model.

- ◆ Weight: 108g(B) / 91g(A)
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0,50mm Pb Strips
- ◆ Frame Material: Nylon
- ◆ Frame Szie: 64-18-150(B) / 69-17-145(A)
- ◆ Color: Black
- ♦ CE: Yes

SK-90 (90 FITOVER)



Designed to fit over prescription eyewear, model 90 fitovers are an excellent choice for those requiring protection against x-ray and liquid splashes.

The large, clear side shields protect the user from lateral exposure while eliminating the "tunnel" effects of look-alikes.

- ◆ Weight: 112g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Standard
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 66-14-120
- ♦ Color: Matte Black
- ♦ C F: Yes

SK-70S (70S ASTRO FLEX) (Without Glasses)





A panoramic goggle constructed of impact resistant propionate.

The frame front has molded side shielding and is fitted to the brow for increased splash protection.

The temple length and frame angle are adjustable for a uniform custom fit

Model 70S comes with Lead side shields for necessary protection from oblique angle exposure.

- Weight: 88g
- ◆ Front Lens: 0.75mm Pb Glass
- Side Shields: 0.50mm Pb Strips
- Frame Material: Propionate
- Frame Szie: 54-19-Adj
- ◆ Color: Clear/Blue
- ◆ CE: Yes

Sung Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmedite	
27		Tel: 82-31-855-7548	Fax: 82-31-855-7549

SK-RE-9941A (9941 ULTRALITE)



The 9941A Ultralites come standard with a rubberized nose bridge. They are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields. The 9941A is a popular frame and ideal for those who have a high power prescription.

- ♦ Weight: 64g
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0.50mmPb Strips Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Szie: 55-22-130
- ◆ Color: Silver
- ♦ CE: Yes

SK-N-BR (NIKE BRAZEN)



The Nike Brazen features an advanced ventilated nose bridge for grip, comfort, and reduced fogging.

Its lightweight nylon design combine streamlined style and enhanced comfort.

The premium Nike Brazen frame offers versatility with its superior coverage and protection.

- ♦ Weight: 75g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 60-16-140
- ◆ Color: Matte Black
- ♦ CE: Yes

SK-RE-RAZER (RAZER)





The Razer's blends a sleek design into a lightweight and durable frame. Its premium features such as a ventilated nose bridge, rubberized temples, and wide-angle lenses maximize comfort and protection from x-ray.

The modern style and use of light-weight materials make the Razer an excellent choice for longprocedures.

Features: Ventilated Nose Bridge, Rubberized Temples for Grip and Comfort, Modern Matte Color Finishes, High Rx Options

- ♦ Weight: 62.5g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb attached
- ◆ Frame Material: Polycarbonate
- ◆ Frame Size: 58-17-135
- ◆ Color: Matte Green, Matte Gray
- ♦ CE: Yes

 Sung Kwang Meditech Co.,Ltd
 www.skmeditech.com
 Contact to E-mail: skmeditech@skmeditech.com

 Tel: 82-31-855-7548
 Fax: 82-31-855-7549
 22

Import Non Lead Aprons

Lead Free

Our lighter lead-free material (0.25 / 0.35 / 0.5) • Environmentally Safe

VEST-SKIRT COMBINATION





This Vest-Skirt combination is the preferred choice for complete protection and weight distribution.

The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue.

Skirt features no-slip, quick release belt for added support.

Double padded shoulders for maximum comfort.

Chest pocket standard.

- (Pictured with Thyroid Tether.)
- Protection : .50mm front / .25mm back
- ◆ Closures : Velcro and Quick-Release Buckle

VEST-SKIRT COMBINATION

APVS

WRAP AROUND BACK RELIEF



Offering full protection, this Back Relief design transfers the weight from the shoulders to the hips with the aid of a 6" bonded elastic belt. (See features of our Back Relief belt.)

This apron features Velcro and buckle closures, as well as a Velcro Breakaway side seam for freedom of movement.

Double padded shoulders for added comfort and a chest pocket are standard features.

- Protection : .50mm front / .25mm back
- ◆ Closures : Velcro and Quick-Release Buckle

WRAP AROUND BACK RELIEF

APWABR

BACK RELIEF





Designed to transfer the weight from the shoulders to the hips.

The apron features a unique padded belt system constructed with a bonded elasticized panel.

The stretch proof panel will not lose its shape while providing increased comfort and back support through the life of the apron.

Extra thick shoulder pads, chest pocket and shoulder strap are standard. Optional tether on shoulder strap available upon request.

- ◆ Protection : .50mm front
- ◆ Closures : Velcro and Quick Release Buckle

BACK RELIEF APBR

APBR

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

FLEXBACK





VEST-SKIRT FLEXBACK





BUCKLE APRON





FLEXBACK w/SURGICAL DROP SHOULDER





Our Flexback apron is one of the most comfortable models you'll ever wear!

Criss-crossing back panels lift the weight off the shoulders and support the lower back.

Optional shoulder strap available upon request.

Extra thick shoulder pads are standard.

Chest pocket standard.

FLEXBACK

Shown with Thyroid Collar.

◆ Protection : .50mm front ◆ Closures : Velcro

APFB

This design combines the benefits of a 2-piece apron with the support of our elasticized flexback apron.

You control the amount of back support by adjusting the Inner elasticized quick release belt.

Chest pocket standard. Shown with Page Cap.

- Protection : .50mm front / .25mm back
- ◆ Closures : Velcro and Quick Release Buckle

VEST-SKIRT FLEXBACK

APVSFB

APBA

Quick Release Buckle Apron is user friendly with large support straps criss-crossing the back allowing proper placement.

Belt is adjustable in order to accommodate varying sizes. Chest pocket and shoulder pads are standard.

Optional shoulder strap (as shown on Back Relief APBR above) available upon request.

- ◆ Protection : .50mm front
- ◆ Closures : Quick Release Buckle

BUCKLE APRON

Two great aprons in one! This is the perfect choice for those who want the convenience of a Surgical Drop Apron and the comfort of our Flexback Apron.

Crisscrossing back panels transfer the weight from the shoulders and distribute it to the hips.

Lift the Velcro tabs for quick removal!

Extra thick shoulder pads and chest pocket are standard. Shown with Apron Cap.

Protection : .50mm front	♦ Closures : Velcro
--------------------------	---------------------

FLEXBACK w/SURGICAL DROP SHOULDER

APFBSD

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Import Non Lead Aprons

SURGICAL DROP-OFF





TIE APRON



VELCRO ADJUSTABLE





MATERNITY APRON





Fluoro done Remove your apron quickly and easily without breaking the sterile field.

Extra thick shoulder pads for added comfort. Chest pocket and hanging straps are standard.

- ◆ Protection : .50mm front
- ◆ Closures : Velcro

SURGICAL DROP-OFF

Our most basic apron. Easy to put on. Ties are long enough to fit a large spectrum of sizes. Chest pocket standard.

- ◆ Protection : .50mm front
- ◆ Closures : Ties

TIE APRON

Crisscrossing back straps keep the shoulders from slipping off. Velcro tabs make adjustability quick and easy.

Extra thick shoulder pads, for comfort and chest pocket are standard.

- ◆ Protection : .50mm front
- ◆ Closures : Velcro

VELCRO ADJUSTABLE

APVA

APSD

APEA

This apron features an adjustable Quick Release Buckle to secure the back panel and easy ties for the front panel allowing complete adjustability during pregnancy. Velcro release shoulders makes removal easy. Apron also includes an additional .50mm equivalent panel.

You can select either the Velcro attachment to the marsupial Inner pouch or Velcro belt.

- Protection : .50mm front / 1.0mm(with added panel) and .25 back
- ◆ Closures : Velcro

MATERNITY APRON

APMA

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

SK-16-1 (Size: 6~9)

GUARD - RR2, RR3



Our Soft-Touch finish creates a smooth interior reducing skin irritation. Anatomically correct curved fingers reduce hand fatigue.

Sizing range from 6 to 9 in half sizes for an unequaled comfortable fit. Attenuates up to 50% more direct beam radiation than other leading products. ThInner fingertips provide enhanced tactile sensitivity.

ATTENUATION

RR-2	RR-2:0.012"/0.30 mm			m RR-3 : 0.024" / 0.60 mm				mm
Kvp	60	90	120		Kvp	60	90	120
Cuff	56%	44%	37%		Cuff	76%	63%	55%
Palm	57%	45%	38%		Palm	75%	62%	54%
Finger	57%	56%	38%		Finger	75%	62%	54%

SKN-16-1 (Size: 6~9)



ATTENUATION

MODE	MODEL 1 : 0.009" / 0.22 mm			MODEL 2 : 0.012" / 0.30 mm				
Kvp	60	90	120	Kvp	60	90	120	
Cuff	51%	40%	34%	Cuff	62%	51%	43%	
Palm	55%	44%	37%	Palm	63%	51%	43%	
Finger	55%	44%	37%	Finger	61%	49%	42%	

SK-SL (Size: S/M/L)



XGUARD Radiation Reducing Sleeve, when used in conjunction with the XGUARD Radiation Reducing Procedure Glove creates an attenuative barrier from your fingertips to your biceps from scatter radiation.

- Inner Nitrile coating reducing proteins, antigens and irritation.
- Lightweight material reduces arm fatigue.
- Extreme flexibility does not restrict movement or dexterity.
- Available in three sizes for a comfortable fit.

ATTENUATION							
XGUARD	60 Kvp	80 Kvp	100 Kvp	120 Kvp	THICKNESS		
Sleeve	55%	45%	35%	31%	0.012" - 0.30mm		

BREAST SHIELD



Sung Kwang Meditech Co.,Ltd

SURGICAL DROP **SHOULDER**



THYROID COLLARS MODEL: AP-TCO




BREAST SHIELD	APSTOLE
SURGICAL DROP SHOULDER	APDRSH
THYROID COLLERS	APTCO / APTCS


www.skmeditech.com

Import Non Lead Aprons Color Swatches




Sung Kwang Meditech $C \in$







Sung Kwang Meditech Co.,Ltd





T) 82-31-855-7548 M) 82-10-4590-7548 F) 82-31-855-7549 E) skmeditech@skmeditech.com Address) 1014-20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea Sung Kwang Meditech co.,Itd

www.skmeditech.com