


Sung Kwang Meditech

**For surgical patient
& Rehabilitation
Radiology Products**


Operating Room

Snow pac (Bean Bag)
Several device
Position pad
Gel product
Instrument Case
Silicon product




Rehabilitation

Gel wheelchair cushion pad
Gel wheelchair back cushion pad



Silver

Gel embo cushion mattress



Radiation Protection For:

Cardiac Cath Lab
Interventional Radiology
Pain Management



Operating Room


SNOW PAC (Bean Bag)

FOR SURGERY POSITIONING PATIENTS

Snow Pac are soft pads (Available in various sizes and shapes) that are filled with thousands of tiny plastic beads. Snow Pac provides more stable and comfortable position to patients and saves valuable time for the entire O.R team on every case.

It prevents pressure points (With Gel). It aids faster, easier set-up and repositioning during surgery.
It aids to become comfortable position. It aids better surgery position presentation.

PRINCIPLE OF SYSTEM



Before

Fix or put the Snow Pac to the patient's body according to the position you want to form (e.g., Supine Position) and connect the suction line to the entrance of the Snow Pac.

Depending on the intensity of the Snow Pac you want to make, use a suction line to inhale the air in the Snow Pac properly. (If you inhale for about 10 to 15 seconds, vacuum is created inside the Snow Pac and it becomes hardens)



After

The inside of the Snow Pac is vacuumed and fixed to the desired shape.

When the suction line is disconnected, the Snow-Pac retains its firm, molded shape and holds the patient securely in position.

After surgery, if user open the entrance of the Snow Pac and press it slightly, air re-enters into the Snow Pac and it becomes soft again and re-useable.

EFFECT

Stable

Because Snow-Pacs are molded well, firmly and securely to shape of the patient's body than sandbags, towel rolls and other devices that are merely propped against the patient.

Comfortable

The apparent wrinkles are inverted-drawn inward by suction and will not irritate the skin.

Better Presentation

With Snow-Pacs the patient can be positioned exactly as the surgeon desires for excellent presentation of the operative site. Snow-Pacs conform to any shape, so the patient's position does not have to be compromised.

More Versatile

With Snow Pacs patient can often be placed in positions that could not be easily achieved using other positioning devices.

No Pressure Points (With Gels)

A Snow-Pac supports all parts of the patient's body uniformly when used with gel like a floating mattress.

Saves Time

Patients can be positioned much more quickly and easily using Snow Pacs, and if necessary, be repositioned during a case.

OTHER FEATURES

Any Vacuum Source

Use wall suction, aspirator pump or any other convenient O.R. vacuum source.

Reusable

As the exterior material is soft and sturdy, Snow Pac can be used for many years as long as they are thoroughly controlled for contamination.

Radiolucent

Excellent radio transparency with only minimal and clearly identifiable artifacts.

Clean, Disinfect, Sterilize

Can be cleaned with regular soaps and liquid disinfectants or sterilized by gas.

Saves Money

Snow Pacs save the cost of sandbags, towel rolls and other positioning devices, and they save time in every case.

Repair Service

Factory inspection and repair service is available at nominal cost after warranty period.
Full One Year Warranty.

INDICATION

KIDNEY POSITION, LATERAL POSITION



A type



A-1 type
A-1-1 type



A-1-2 type
A-1-3 type

Easily up in nephrectomy (A type), and also for handicapped person with cerebral palsy.(A-1 type).



A type	73(L) * 100(W)cm
A-1 type	73(L) * 100(W)cm
A-1-1 type	90(L) * 100(W)cm
A-1-2 type	73(L) * 100(W)cm(U)
A-1-3 type	90(L) * 100(W)cm(U)

THYROID POSITION



B type
B-1 type

Provides hyperextension, prevents postoperative discomfort.

B type	72(L) * 88(W)cm
B-1 type	118(L) * 88(W)cm

HEAD STABILIZATION POSITION



C type

Stabilizes head in eye, facial or dental surgery.

C type	30(L) * 40(W)cm
--------	-----------------

SUPINE POSITION



D type




Supports and eliminates back pains during long procedures.

Prevents foot drop splint in Emergency.

D type	30(L) * 55(W)cm
--------	-----------------

VARIOUS FULL BODY POSITION



E type

E type	165(L) * 72(W)cm
--------	------------------

POSITION PADS

PAD FOR SUPINE POSITION




Special sponge (Inner)
Good positioner
Anti-bedsores effect Very soft, Good elasticity Good endure.
Poly urethane vinyl (Outer)
Water resistant.

SUPINE

Head, foot pad	
SK-1	22(L) * 16(W) * 11 (5.5) (H)cm
Head pad	
SK-1-1	22.8(O.D) * 11.5(I.D) * 5(H)cm
Shoulder pad	
SK-2	23(L) * 14(W) * 7(H)cm


PAD FOR PRONE POSITION



PRONE

Chest pad (thick, thin)	
SK-5-1	18(L) * 52(W) * 13(H)cm
Head pad	
SK-6	32(L) * 25(W) * 10(H)cm
SK-6-1	24(L) * 26.5(W) * 18(H)cm
Knee pad	
SK-7	30(L) * 20(W) * 6(H)cm
Leg pad	
SK-3	40(L) * 17(W) * 17(H)cm

PAD FOR LATERAL POSITION



LATERAL

Arm pad	
SK-4	12(L) * 62(D) * 7(H)cm
Leg pad	
SK-8	27(L) * 10(D)cm
SK-8-1	45(L) * 9(D)cm
SK-8-2	45(L) * 11 (D)cm
SK-8-3	47(L) * 13.5 (D)cm
SK-9	50(L) * 25(W) * 7(H)cm
SK-10	65(L) * 40(W) * 6(H)cm
Shoulder pad	
SK-12	22(L) * 15(W) * 12(H)cm


Autoclave Instrument Cases

Instrument Sterilization Cases are designed to protect fine surgical instruments and various scopes during sterilization and storage. Our cases are made of elastic special materials for flexible handling and **it is transparent so that the inside can be seen during autoclaving**. There are fine circular holes in the upper(Lid) and lower(Base) of the case to **help circulate steam during sterilization**. Silicon Embo Mat or Bar is placed on the bottom and prevents instruments from moving freely.

- This sterilizing instrument tray is designed to withstand repeated exposure to **134 degrees (Steam / EO / Plasma sterilizable)**.
- There are a variety of designs and sizes suitable for device applications such as **micro instrument, Laparoscope, Telescope, and Endoscope**.
- The size table below is based on the outer diameter, and the inner diameter by -1 to 2cm.


Micro Instrument Case(26(L) x 17(W))

SK-19(2F)	26(L) x 17(W) x 4(H) cm
SK-19-a(1F)	26(L) x 17(W) x 2.5(H) cm



Scope Case(Telescope, Laparoscope)

SK-20-1 (Laparoscope Case)	48(L) x 9(W) x 5(H) cm
SK-20-2(Telescope Case)	28(L) x 7(W) x 5(H) cm



Micro Instrument Case(1F)

SK-19-S	16(L) x 6(W) x 3(H) cm
SK-19-1	25(L) x 5(W) x 2.5(H) cm
SK-19-3	30(L) x 10(W) x 2.5(H) cm



Micro Instrument Case(1F, 2F)

SK-19-5(2F)	39(L) x 27(W) x 4(H) cm
SK-19-5-a(1F)	38(L) x 26(W) x 2.5(H) cm
SK-19-5-b(1F)	39(L) x 27(W) x 12(H) cm




Scope Case(Long, Big size)

SK-20-3(Endoscope)	37(L) x 15(W) x 7(H) cm
SK-20-4(Laparoscopic)	60(L) x 20(W) x 6(H) cm
SK-20-5(Long scope)	72(L) x 15(W) x 10(H) cm
SK-21-1 (Laparoscopic)	64(L) x 25(W) x 15.5(H) cm



Ophthalmic Instrument Case(Silicon color inside: blue)

SK-19-a(1F)	26(L) x 17(W) x 2.5(H) cm
SK-19-S	16(L) x 6(W) x 3(H) cm
SK-19-5-a(1F)	38(L) x 26(W) x 2.5(H) cm




CASES

INSTRUMENT STERILIZATION TRAYS USE CASES

The instrument sterilization case photos below are examples of the use of cases that were actually applied in storage or sterilization processes by various departments such as surgery and internal medicine.




SK-19-5(2F)
39(L) x 27(W) x 4(H) cm



SK-20-3
37(L) x 15(W) x 7(H) cm



SK-19-a(1F)
26(L) x 17(W) x 2.5(H) cm



SK-19-5-b(1F)
39(L) x 27(W) x 12(H) cm




SK-21-1
64(L) x 25(W) x 15.5(H) cm



SK-19-3
30(L) x 10(W) x 2.5(H) cm



SK-20-3
37(L) x 15(W) x 7(H) cm



SK-19-1
25(L) x 5(W) x 2.5(H) cm



SK-20-2
28(L) x 27(W) x 5(H) cm



SK-19-5-a(1F)
38(L) x 26(W) x 2.5(H) cm

Strengths of Sungkwang Meditech trays:

- Cases are transparent. (Users can see instruments during autoclaving)
- Cases are light but strong.
(We do not use stainless, steel materials. But, use special strong materials)
- There are fine circular holes in the upper(Lid) and lower(Base).
(It helps circulate steam during sterilization)
- Silicon Embo Mat or Bar is available.
(It is placed on the bottom and prevents instruments from moving freely).
- It is possible to assemble freely based on the shape of the instrument.




SK-21-1
64(L) x 25(W) x 15.5(H) cm




SK-20-3
37(L) x 15(W) x 7(H) cm



SK-20-5
72(L) x 15(W) x 10(H) cm




SK-20-4
60(L) x 20(W) x 6(H) cm



SK-20-4
60(L) x 20(W) x 6(H) cm




SK-20-3
37(L) x 15(W) x 7(H) cm



SK-20-3
37(L) x 15(W) x 7(H) cm



SK-20-3
37(L) x 15(W) x 7(H) cm



SK-20-4
60(L) x 20(W) x 6(H) cm

Operating Room

CASES

NEEDLE COUNTER & SHARPS COLLECTOR

Inner : magnet
For surgery in OR
It aids to count needle thoroughly.




SK-18



SK-18-1

NEEDLE HUNTER

It finds needles by rolling on the floor



SK-13-2

(External)

SK-18 (S)	14(L) * 4(W) * 1.5(H)cm
SK-18-1 (L)	17(L) * 8(W) * 2(H)cm
SK-13-2	100(L) * 30(W)

MICRO INSTRUMENT ONE FLOOR CASE



SK-19-a (1F)



(External)

SK-19-a (1F)	26(L) * 17(W) * 2.5(H)cm
--------------	--------------------------

MICRO INSTRUMENT TWO FLOOR CASE



SK-19 (2F)



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)

SK-19(2F)	26(L) * 17(W) * 4(H)cm
-----------	------------------------

MICRO INSTRUMENT ONE FLOOR CASE



SK-19-S

For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.



SK-19-1

Steam / EO / Plasma Sterilizable




SK-19-3

(External)

SK-19-S	16(L) * 6(W) * 3(H)cm
SK-19-1	25(L) * 5(W) * 2.5(H)cm
SK-19-3	30(L) * 10(W) * 2.5(H)cm

Operating Room

MICRO INSTRUMENT ONE FLOOR CASE



SK-19-5-a



SK-19-5-b




Steam / EO / Plasma Sterilizable

(External)

SK-19-5-a(1F) 38(L) * 26(W) * 2.5(H)cm

SK-19-5-b(1F) 39(L) * 27(W) * 12(H)cm

MICRO INSTRUMENT TWO FLOOR CASE



SK-19-5 (2F)

For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.


Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)

SK-19-5(2F) 39(L) * 27(W) * 4(H)cm

LAPAROSCOPE CASE



For Laparoscope.

Silicon bar is placed to prevent laparoscope to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)	
SK-20-1	48(L) * 9(W) * 5(H)cm

TELESCOPE CASE



For Telescope.

For Arthroscope Silicon bar is placed to prevent laparoscope to move freely.


Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)	
SK-20-2	28(L) * 7(W) * 5(H)cm

Operating Room

ENDOSCOPE INSTRUMENT CASE



SK-20-3

For Endoscope Instrument.

Silicon bar is placed to prevent laparoscopic instrument to move freely.

Cases are designed to with stand repeated exposure to 134°C steam temperature.

SK-20-3 is used for endoscope camera case.

Steam / EO / Plasma Sterilizable

(External)
SK-20-3 37(L) * 15(W) * 7(H)cm

LAPAROSCOPIC INSTRUMENT CASE



SK-20-4

Silicon bar is placed to prevent laparoscopic instrument to move freely.

Cases are designed to with stand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable


It contains 2pcs or more of laparoscopic instrument.

(External)
SK-20-4 60(L) * 20(W) * 6(H)cm

LONG SCOPE INSTRUMENT CASE



SK-20-5



Sterilization cases used during laparoscopic long scope and robotic instrument surgery.

Silicon bar or mat is placed to prevent movement of instruments used in obstetrics, urology, otorhinolaryngology, surgery and robotic surgery.

Two or more laparoscopic instruments can be inserted.



Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)

SK-20-5 72(L) * 15(W) * 10(H)cm

LAPAROSCOPIC INSTRUMENT CASE



SK-21-1

Silicon bar is placed to prevent laparoscopic instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable


It contains 10pcs of laparoscopic instrument.

(External)

SK-21-1 64(L) * 25(W) * 15.5(H)cm

Operating Room

OPHTHALMIC INSTRUMENT CASE



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable


(External)

SK-19-S	16(L) * 6(W) * 3(H)cm
SK-19-a (1F)	26(L) * 17(W) * 2.5(H)cm
SK-19-5-a (1F)	38(L) * 26(W) * 2.5(H)cm

BIOPSY CASE

R2	Delphain	Pyramidal	L2
R3			L3
R. Jungular	R. Thyroid	L. Thyroid	L. Jungular
R4			L4
R5	R.CN	L.CN	L5
R. Para. T			L. Para. T

SK-21-a-1 (for Thyroid)



Material : stainless steel with marked L/N number,

Use : for specimen storage

order made available

Rt Level 1		Rt Level 1	
Rt Level 2		Rt Level 2	
Rt Level 3		Rt Level 2	
Rt deep margin	Rt rotter's	Lt deep margin	Lt rotter's
Rt. axillary	Rt. sentinel	Lt. axillary	Lt. sentinel

SK-21-a-2 (for breast)



(External)

SK-21-a-1	for Thyroid 16.5(W) * 19(L) * 3.2(H)
SK-21-a-2	for breast 14.5(W) * 15(L) * 3.2(H)
SK-21-a-3	for stomach 31.5(W) * 8.5(L) * 3.2(H)
SK-21-a-4	for biliary and Pancreas 30(W) * 8.5(L) * 3.2(H)


1	2	3	4	5	6	7	4sa		
8	9	10	11	11p	11d	12	14v	4sb	4d

SK-21-a-3 (for stomach)



5	6	8	12a	12b	12c	12h	12p
13	14a	14v	16	16a2	16b1	17	

SK-21-a-4 (for biliary and Pancreas)



GEL PRODUCTS

Capillary closure is at the root of ulcer formation.

The blockage or damage of capillaries prevents the natural maintenance of healthy skin.


In operation, the use of proper positioning devices prevent excessive pressure build-up and capillary closure.

Evenly distribute weight and provide constant support beneath patients.

Gel products are cost effective, latex free, fast & easy cleaning.



SK-23



SK-24



SK-24-1



SK-25



SK-29



Pressure Injury Free Gel Products



SK-22-1



SK-26



SK-22-2



SK-22




SK-29



SK-28



SK-27



DOUGHNUT GEL HEAD PAD



Doughnut gel head pad provides excellent protection when the patient is positioned for surgery by supine position.

SK-22(Adult)	20(D) * 4.5(H)cm
SK-22-3(Child)	15(D) * 3.5(H)cm
SK-22-4(Baby)	9.2(D) * 2.4(H)cm

C-DOUGHNUT GEL HEAD PAD



C-doughnut gel head pad provides clear air passage way when the patient is positioned for surgery by lateral or prone position.

SK-22-1 (Adult)	21 (D) * 4.5(H)cm
SK-22-1-1 (Child)	15(D) * 3.5(H)cm
SK-22-1-2 (Baby)	9.2(D) * 2.4(H)cm

GEL OPHTHALMIC HEADREST



It was designed for stable head placement for all ophthalmic and oral or facial procedures.

Also provides patient with comfort while under local anesthesia, minimizing patient movement.

It was also designed for surgeon's functional arm support.

SK-22-a(Adult)	31 (L) * 27 (W) * 10(H)cm
SK-22-b(Child)	31 (L) * 25 (W) * 7.3(H)cm

GEL OPHTHALMIC CRADLE HEADREST




It was designed for ENT, plastic surgery, all ophthalmic cases.

It helps prevent risk of alopecia especially during long surgery.

Cradling of patient's head provides stability.

SK-22-c(Adult)	23 (L) * 23 (W) * 5.7(H)cm
SK-22-d(Child)	19 (L) * 21.6 (W) * 3.8(H)cm

GEL HEAD PAD



It provides optimum stability of patient's head while in the supine position.

Can also be used for position and protect the knee or calf.

SK-22-e(Adult)	18 (L) * 18 (W) * 8(H)cm
SK-22-5(Baby)	16 (L) * 12 (W) * 1(H)cm

PRONE GEL HEAD PAD



It provides clear air passage way when the patient is positioned for surgery by prone position.

SK-22-2	27 (L) * 24 (W) * 15/12(H)cm
SK-22-2-1	27 (L) * 21 (W) * 8(H)cm


GEL BODY PAD



Gel body pad provides complete pressure reduction and shares across entire table surface

SK-23-2	84(L) * 50(W) * 1.3(H)cm
SK-23	160(L) * 50(W) * 1.3(H)cm
SK-23-1	183(L) * 50(W) * 1.3(H)cm

GEL BODY PAD 3-PIECE SET



This set includes head, torso and foot section pads for use with segment tables.

SK-23-1-1 Gel body pad 3-piece set(a+b+c)	
SK-23-1-1a(Head)	50(L) * 23(W) * 1.3(H)cm
SK-23-1-1b(torso)	100(L) * 50(W) * 1.3(H)cm
SK-23-1-1c(foot)	50(L) * 50(W) * 1.3(H)cm


GEL INCUBATOR PAD



It provides gentle pressure for neo-natal patients in the NIC units. It is designed for the isolette beds.

SK-23-5	60(L) * 33(W) * 1.3(H)cm
SK-23-5-1	60(L) * 40(W) * 1.3(H)cm


GEL CHAIR PAD



Gel chair pad is used underneath prominence area in a variety of applications.

SK-24	42(L) * 42(W) * 1.3(H)cm
SK-24-2	38(L) * 38(W) * 1.3(H)cm
SK-24-2-1	40(L) * 40(W) * 1(H)cm
SK-24-2-2	40(L) * 45(W) * 1(H)cm
SK-24-3	32(L) * 32(W) * 1.3(H)cm

GEL EMBO PAD




Gel embo chair pad is used for iliac crest, under heels and elbows. Is used as kidney rest, axillary roll while patient in lateral position.

Gel embo body pad provides complete pressure reduction and shears cross entiret ablesur face.

Gel embo chair pad	
SK-24-1	40(L) * 45(W) * 1.6(H)cm
Gel embo body pad	
SK-23-3	160(L) * 45(W) * 1.6(H)cm
Gel embo half body pad	
SK-23-4	80(L) * 45(W) * 1.6(H)cm

Operating Room

GEL KNEE PAD



SK-25 (S)
SK-25-1 (L)

Gel knee pad is used for OBGY surgery in lithotomy position. It provides against nerve injury.


SK-25	40 (L) * 24 (W) * 1 (H)cm
SK-25-1	50 (L) * 52 (W) * 1 (H)cm

GEL LATERAL PAD GEL



SK-26

GEL THYROID PAD



SK-26-2

Gel lateral pad is designed to accommodate both the shoulder and provides exceptional pressure relief and nerve protection.

Gel thyroid pad is designed for thyroidectomy. It provides patient with gentle pressure while long surgery with hyperextension position. The tissue friendly feels make it very comfortable.

SK-26	70 (L) * 50 (W) * 12 (H) / 0.8cm
SK-26-2	45 (L) * 30 (W) * 7 (H)cm

GEL PATIENT POSITIONER



SK-27-a

SK-27-b

SK-27-c

It is good for any procedure. versatile positioner.

SK-27-a (Small)	25 (L) * 13 (W) * 3 (H)cm
SK-27-b (Large)	35 (L) * 12 (W) * 3 (H)cm
SK-27-c (Extra long)	46 (L) * 11 (W) * 3 (H)cm

GEL PRONE CHEST PAD



SK-27-d

SK-27-e

They were designed for uniformly conform to a patient's shoulder to chest area, with rounded edges for the axilla.

They provide optimum alleviating pressure on the thorax area.

SK-27-d (Right)	70/58 (L) * 20 (W) * 5 (H)cm
SK-27-e (Left)	70/58 (L) * 20 (W) * 5 (H)cm

GEL SHOULDER PAD (CHEST)



SK-27



SK-27-1 (hard)



SK-27-3



SK-27-2 (hard)



SK-27-4

Semi-circular design provides maximum pressure protection while patient is in supine position.

SK-27, 27-4 can be used under the knee, the foot against the board.

While patient is in lateral position, SK-27 can be used for the shoulder protection.

In thyroidectomy, SK-27-4 is used for neck hyperextension.

while patient is in prone position SK-27-1, SK-27-2 can be used under the chest.

SK-27	47 (L) * 10 (W) * 5 (H)cm
SK-27-1	47 (L) * 15 (W) * 8 (H)cm
SK-27-2	47 (L) * 15 (W) * 12 (H)cm
SK-27-3	25 (L) * 15 (W) * 8 (H)cm
SK-27-4	52 (L) * 10 (W) * 8 (H)cm

GEL HEEL & WRIST PAD



SK-28

Gel heel & wrist pad provides safety and protects heel & wrist area body prominence.

SK-28	20(L) * 11(W) * 7(H)cm
-------	------------------------

GEL HEEL SUPPORTER



SK-28-1

GEL HEEL CUP



SK-28-2

It protects heel area without hook and loop attachment.

It protects the patient's Achilles tendon area.

It supports heels. Appropriate for recovery.


SK-28-1	16(L) * 8.5(W) * 5(H)cm
---------	-------------------------

SK-28-2	17(L) * 18(W) * 1.3(H)cm
---------	--------------------------

GEL HEEL PROTECTOR



(cover)



SK-28-3

It supports heels appropriately for recovery.

SK-28-3	17(L) * 23(W) * 1.3(H)cm
---------	--------------------------

GEL ARMBOARD PAD



SK-29



It is applied for supine position for protection of bony prominence area.

SK-29	55(L) * 15(W) * 1.3(H)cm
-------	--------------------------

GEL LATERAL ARMBOARD PAD



SK-29-1



SK-29-1a

It is applied for lateral position.

Designed to provide arm support and stability.

SK-29-1 (Small)	41(L) * 15(W) * 5(H)cm
-----------------	------------------------

SK-29-1a (Large)	55(L) * 15(W) * 6(H)cm
------------------	------------------------

GEL WRIST AND ANKLE PROTECTOR GEL ARMBOARD BINDER



SK-17-a (S)



SK-17-b (L)



SK-17-c

It is used for nerve protection in fastening wrist and ankle.

SK-17-a Gel Wrist protector	26(L) * 6(W) * 1(H)cm
-----------------------------	-----------------------

SK-17-b Gel Ankle protector	33(L) * 6(W) * 1(H)cm
-----------------------------	-----------------------


SK-17-c Gel Armboard binder	21(L) * 6(W) * 1(H)cm
-----------------------------	-----------------------

GEL WILSON FRAME COVER (2 PIECE SET)

SK-17-1



GEL SK 5-1 COVER (2 PIECE SET)



SK-5-1a

It provides comfort and distributes pressure in prone position spine surgery.

SK-17-1	75(L) * 24.5(W) * 1.3(H)cm
---------	----------------------------

SK-5-1a	52(L) * 31(W) * 1.3(H)cm
---------	--------------------------

GEL RELTON FRAME PAD (4 PIECE SET)

SK-17-1-1



It provides comfort and safety in prone position spine surgery.

SK-17-1-1	37(L) * 17(W) * 14(H) / 1.3(thick)cm
-----------	--------------------------------------

Operating Room

GEL FOREARM PROTECTOR



SK-17-2 (L)
SK-17-2-1 (S)



It is used for nerve protection of arm and elbow in supine position surgery.

SK-17-2 Gel forearm protector(L)

46(L) * 21 (W) * 0.8(H)cm

SK-17-2-1 Gel forearm protector(S)

38(L) * 11.5(W) * 0.8(H)cm

GEL BODY BELT



SK-17-3



It helps safety, stability and tissue protection during surgery.

SK-17-3

50(L) * 8.5(W) * 0.8(H)cm

GEL SHOULDER PROTECTOR



SK-17-3-1

It is designed for patient in trendelenburg down position to protect shoulder areas.

SK-17-3-1

16(L) * 15(W) * 2(H)cm

GEL ALLEN PAD

SK-17-4-1



It is used lithotomy position for nerve protection.

SK-17-4-1

68(L) * 38(W) * 1(H)cm

GEL ANKLE HOLDER (2 PIECE SET)

SK-17-5



It is used in lithotomy position for nerve protection.

SK-17-5

40(L) * 3(W) * 1(H)cm

GEL SNOW PAC A-1 type COVER



SK-17-6



Designed for A-1 type snow pac. It provides comfort.

SK-17-6

100(L) * 73(W) * 1(H)cm

GEL SNOW PAC B type COVER



SK-17-6-1



Designed for B type snow pac.

SK-17-6-1

88(L) * 72(W) * 1(H)cm

GEL SNOW PAC C type COVER



SK-17-6-2



Designed for C type snow pac.

SK-17-6-2


30(L) * 40(W) * 1(H)cm

SILICON

SILICON ROD & MAT

Silicon rod & mat protect microinstruments in microsurgery.


It endures 300 degree C temperature up.



Hole mat



Embo mat



SILICON ROD

SK-31 (Rod)	1.5(L) * 10(W) * 0.8(H)cm
SK-31-1 (Rod)	1.5(L) * 20(W) * 0.8(H)cm
SK-31-1-1 (Rod)	1.5(L) * 20(W) * 2.8(H)cm
SK-31-2 (Rod)	1.5(L) * 30(W) * 0.8(H)cm
SK-31-2-1 (Rod)	1.5(L) * 30(W) * 2.8(H)cm

SILICON MAT

SK-31-3(Mat)	50(L) * 33(W) * 0.2(H)cm
--------------	--------------------------

SILICON HOLE MAT

SK-31-3-a	50(L) * 25(W) * 0.2(H)cm
SK-31-3-b	33(L) * 25(W) * 0.2(H)cm
SK-31-3-c	25(L) * 25(W) * 0.2(H)cm
SK-31-3-d	50(L) * 33(W) * 0.2(H)cm


SILICON EMBO MAT

SK-31-4	23.5(L) * 15(W) * 1.5(H)cm
SK-31-4-1	37(L) * 24.3(W) * 1.5(H)cm

SILIMAGNET MAT

SK-31-5	37.5(L) * 29(W) * 0.6(H)cm
SK-31-5-1	43(L) * 35(W) * 0.6(H)cm

Silimagnet mat



SILICON TUBE

SILI GAS LINE (SK-31-4-3)

In laparoscopic surgery used for CO₂ gas infusion.



SILI SUCTION TUBE (SK-31-4-4)

In neurosurgery, used for suction tube. It provides elasticity, flexibility



SK-31-4-3	5.0(I.D) * 7.0(O.D) mm
SK-31-4-4	6.0(I.D) * 11.0(O.D)mm

Several Devices

PATIENT TRANSPORT BOARD



The fully equipped patent transport board is ideal for transport when patients need to be transferred from MRI, operating table, bed, etc.

SIZE :

SK-MTB80	80 x 35cm
SK-MTB140	140 x 45cm
SK-MTB183	183 x 55cm


RADIATION SHIELDING X-RAY CURTAINS FOR SCANNER



Radiation Shielding Curtains for security, food scanners. The protection range includes lead equivalent values of 0.125, 0.25, 0.375, 0.5mm depending on the level of protection required.

Customer made.

ARM PROTECTOR



SK-13

CUSHION PAD

Multi-Purpose position pillow for hospital operating room.

It is easy to use for short-term surgery like Hemorrhoid surgery

PROTECTIVE BODY BAND

PROTECTIVE KNEE BAND



SK-14



SK-14-1



Inner : soft cotton cushion

SK-CU-CP(S)



SK-CU-CP(L)



WRIST BAND

Elbow Band



SK-14-6



SK-14-7

SK-13	40.5(L) * 17(H) * 15.5(W)cm Color: White
SK-14	200(L) * 11(W)cm
SK-14-1	70(L) * 7(W)cm
SK-CU-CP-L	45(L) * 70(W) cm
SK-CU-CP-S	40(L) * 60(W) cm
SK-14-6	18(L) * 18.5(W) cm
SK-14-7	35(L) * 12(W) cm

Custom made available

GEL WHEELCHAIR CUSHION PAD WITH COVER

It was made with a high-quality polyurethan gel and PU foam.


The inner gel material is highly resistant to bacteria and bodily fluids, and can be easily cleaned with warm soapy water or a disinfectant. Also, it helps dissipate heat buildup.

Gel Wheelchair cushion pad offers excellent pressure reduction & distribution, great comfort, and positioning stability.


It can be applied for individuals with a history of skin breakdown or those who suffer from pressure sores.

Feel free to enjoy our products below.

GEL WHEELCHAIR CUSHION 2"




GEL WHEELCHAIR CUSHION 3"




SK-1616-2	16(L) * 16(W) * 2(H) Inch
SK-1816-2	18(L) * 16(W) * 2(H) Inch
SK-1818-2	18(L) * 18(W) * 2(H) Inch
SK-1616-3	16(L) * 16(W) * 3(H) Inch
SK-1816-3	18(L) * 16(W) * 3(H) Inch
SK-1818-3	18(L) * 18(W) * 3(H) Inch

GEL WHEELCHAIR CUSHION CURVED 3.5"




GEL WHEELCHAIR CUSHION MORE CURVED 4"




SK-C1616-3.5	16(L) * 16(W) * 3.5(H) Inch
SK-C1816-3.5	18(L) * 16(W) * 3.5(H) Inch
SK-C1818-3.5	18(L) * 18(W) * 3.5(H) Inch
SK-MC1616-4	16(L) * 16(W) * 4(H) Inch
SK-MC1816-4	18(L) * 16(W) * 4(H) Inch
SK-MC1818-4	18(L) * 18(W) * 4(H) Inch

GEL WHEELCHAIR EMBO CUSHION 2"




GEL WHEELCHAIR EMBO CUSHION 3"




SK-E1616-2	16(L) * 16(W) * 2(H) Inch
SK-E1816-2	18(L) * 16(W) * 2(H) Inch
SK-E1818-2	18(L) * 18(W) * 2(H) Inch
SK-E1616-3	16(L) * 16(W) * 3(H) Inch
SK-E1816-3	18(L) * 16(W) * 3(H) Inch
SK-E1818-3	18(L) * 18(W) * 3(H) Inch

GEL WHEELCHAIR EMBO CURVED CUSHION 3.5"



GEL WHEELCHAIR EMBO MORE CURVED CUSHION 4"




SK-EC1616-3.5	16(L) * 16(W) * 3.5(H) Inch
SK-EC1816-3.5	18(L) * 16(W) * 3.5(H) Inch
SK-EC1818-3.5	18(L) * 18(W) * 3.5(H) Inch
SK-EMC1616-4	16(L) * 16(W) * 4(H) Inch
SK-EMC1816-4	18(L) * 16(W) * 4(H) Inch
SK-EMC1818-4	18(L) * 18(W) * 4(H) Inch

GEL WHEELCHAIR BACK CUSHIONS WITH COVER

Polyurethane gel and PU foam helps comfort and durability for long-term.

The embossed gel feature aids comfort for long-lasting back and lumbar support in long-term seating environments.

EMBO GEL BACK CUSHION(A type)



WHEELCHAIR BACK CUSHION COVER


SK-WB-cover




SK-B1417-a	14(L) * 17(W) * 2(H) Inch
SK-B1617-a	16(L) * 17(W) * 2(H) Inch
SK-B1817-a	18(L) * 17(W) * 2(H) Inch
SK-B2017-a	20(L) * 17(W) * 2(H) Inch

SK-B1417-b	14(L) * 17(W) * 2(H) Inch
SK-B1617-b	16(L) * 17(W) * 2(H) Inch
SK-B1817-b	18(L) * 17(W) * 2(H) Inch
SK-B2017-b	20(L) * 17(W) * 2(H) Inch

EMBO GEL BACK CUSHION (B type)



EMBO GEL BACK CUSHION(C type)



SK-B1417-c	14(L) * 17(W) * 2(H) Inch
SK-B1617-c	16(L) * 17(W) * 2(H) Inch
SK-B1817-c	18(L) * 17(W) * 2(H) Inch
SK-B2017-c	20(L) * 17(W) * 2(H) Inch

GEL PRODUCTS-ACCESSORY

EMBO GEL WHEELCHAIR ARMBBOARD PAD



Adhesive tape attached in bottom

RH-3	23(L) * 7(W) * 1.6(H) cm
------	--------------------------

WHEELCHAIR CUSHION COVER



Bottom cover is made by non slip material.

Side cover is made by net shaped material.

Top cover is made by very soft, skinlike material.

SK-W-cover	
------------	--

GEL TOILET COVER / GEL COMMODE SEAT



It is very effective for users with existing skin breakdown.

RH-4(Dual)+cover	39(L) * 10(W) * 1(H)cm
RH-6 (close type)	45(L) * 45(W)cm
RH-6-1 (open type)	45(L) * 46(W)cm

SNOW PAC



AIR SUCTION



HB-102a



RH B type

It is used for lurgeric patient and unpositioned handicapped person. Make position easy and fast. Then, suction equipment is used.

RH B type	83(L) * 108(W)cm
HB-102a	33cm

GEL HORSE S EQUESTRIAN PAD

RECTANGULAR
ROUND



SK-EQ-a Rectangular shape with cover

36(L) * 50(W) * 1 (H) cm

SK-EQ-a cover Rectangular shape cover only

36(L) * 50(W) * 1 (H) cm

SK-EQ-b Round shape with cover

40(L) * 55(W) * 1 (H) cm

SK-EQ-b cover Round shape cover only

40(L) * 55(W) * 1 (H) cm

SILVER PRODUTS

GEL EMBO CUSHION PILLOW / MATTRESS WITH COVER (Memory foam + embo gel)


Characteristics:

C-spine supporting effect is developed by medical basis.

As humanengineering structure, it softly supports head and provides comfortable.



RH-1



RH-2

Outer : 100% cotton velvet with silvercholid coating

Inner : cotton 75% polyester 25%

Materials : slowly recovery
polyurethane foam 100%
+ embo gel pad

Color : silver, ivory

Kind : therapeutic, antibacterial,
deodorizational.

RH-1	50(L) * 33(W)cm
RH-2 (twin)	RH-2-1 (full)
RH-2-3 (queen)	RH-2-4 (king)

Cover Only

RH-2-cover	RH-2-1-cover
RH-2-3-cover	RH-2-4-cover

Lead Aprons

Several Types of Aprons For Radiation Protection



SK-15 (Buckle type)



SK-15-1 (Coat type)



SK-15-4 (X-band type)



SK-15-5 (Buckle X-band type)



SK-15-6 (X-buckle type)

Exposure to radiation is fatal to the medical professional and the human body.

Here are all the radiation protection products you need for protection from radiation exposures.

We manufacture a comprehensive range of **high-quality, lightweight, and cost-effective** products for medical customers around the world.

Inner : Natural rubber-based Lead flexible wearing comfort and best.

Outer : nylon fabric, etc satisfies the workmanship hygienic requirements and the demands of medical use.

Pb Equivalent :
0.25 / 0.35 / 0.50 mmPb Available

±2~4cm, unit : cm
L : Length, W : width

MALE

	S		M		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	93	58	98	60	103	63	108	66
Coat type (Front/back)	96/86	53	99/89	55	105/95	57	106/96	60
X-band type	94	55	96	57	100	60	104	62
Vest/Skirt type	59/50	49/42	60/52	51/43	61/53	53/46	63/54	56/49
panorama type	82	47	85	50	88	53	91	55

FEMALE

	S		M		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.

Available to
SK-15, SK-15-1, SK-15-4, SK-15-5

Several Types of Aprons For Radiation Protection



SK-15-2 (Vest & Skirt type)



Lead Free



SK-15-3-2

Thyroid shield is an essential to protect the thyroid gland from radiation exposure. (SK-16, SK-16-a, SK-16-b)

Also, We provides thyroid shield hygiene covers for cleanliness. (SK-601, 602)

Lead gloves series provide full protection from finger to arm (SK-500, 500-a, 500-b)

Find a wide variety of radiation protection products.

Pb Equivalent :

0.25 / 0.35 / 0.50 mmPb Available



SK-15-3 (Half type)



SK-15-3-1 (Panorama type)



Color type



Blue Cross



Dark Blue



Pink



SK-500



SK-16




SK-16-a



SK-16-b



SK-601
(Reusable cover)



SK-602
(Disposable cover)



SK-500-a



SK-500-b



SK-16-4



SK-16-4-a



SK-16-5



SK-16-5-a

Our Lead cap and gonad shield are **lightweight** and **comfortable**.

An adjustable velcro in the back is attached for a **customizable fit**.

-Two types for cap : SK-16-4 / SK-16-4-a

-Two types for gonad : SK-16-5 / SK-16-5-a

-Order Types : Male, female

-Protection : 0.25mm front

-**Pb Equivalent** : 0.25 / 0.35 / 0.50 mmPb Available

Custom-made available.

Lead cap	SK-16-4 (Basic Model)
Lead cap	SK-16-4-a (Cool wear moisture-wicking fabric)
Lead gonad shield	SK-16-5
Lead gonad shield	SK-16-5-a

Non Lead Aprons

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15 (Buckle type)



SKN-15-1 (Coat type)



SKN-15-4 (X-band type)



SKN-15-5 (Buckle X-band type)



SKN-15-6 (X-buckle type)

Exposure to radiation is fatal to the medical professional and the human body.

Here are all the radiation protection **Lead-Free** products you need for protection from radiation exposures.

We manufacture a comprehensive range of **high-quality, lightweight, and cost-effective** products for medical customers around the world.

Inner : Natural rubber-based Lead flexible wearing comfort and best.

Outer : Water Resistant, Abrasion Resistant, Cool wear

Pb Equivalent : 0.25 / 0.35 / 0.50 mmPb Available

±2~4cm, unit : cm
L : Length, W : width

MALE

	S		M		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	93	58	98	60	103	63	108	66
Coat type (Front/back)	96/86	53	99/89	55	105/95	57	106/96	60
X-band type	94	55	96	57	100	60	104	62
Vest/Skirt type	59/50	49/42	60/52	51/43	61/53	53/46	63/54	56/49
panorama type	82	47	85	50	88	53	91	55

FEMALE

	S		M		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.

Available to
SKN-15, SKN-15-1, SKN-15-4, SKN-15-5

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15-2 (Vest & Skirt type)



SKN-15-3 (Half type)



SKN-15-3-1 (Panorama type)



SKN-15-3-2



SKN-16



SKN-16-a

SKN-16-b

Various types of aprons available.


Thyroid shield is an essential to protect the thyroid gland from radiation exposure.
(SKN-16, SKN-16-a, SKN-16-b)

Find a wide variety of radiation protection products and colors.


Pb Equivalent :

0.25 / 0.35 / 0.50 mmPb Available

Color type




Dia Navy




Dia Orange



Dia Chocolate



Pearl Grey



Pearl Navy



Pearl Purple



Pearl Green



Pearl Gold



Melange Purple



Melange Orange



Melange Navy



Black Dot



Honeycomb Grey



SKN-16-4



SKN-16-4-a



SKN-16-5



SKN-16-5-a

Our Lead-Free cap and gonad shield are **lightweight** and **comfortable**.

An adjustable velcro in the back is attached for a **customizable fit**.

-Two types for cap : SKN-16-4 / 16-4-a

-Two types for gonad : SKN-16-5 / 16-5-a

-Order Types : Male, female

-Protection : 0.25mm front

-**Pb Equivalent** : 0.25 / 0.35 / 0.50 mmPb Available

Custom-made available.

Lead Free cap SKN-16-4 (Basic Model)


Lead Free cap SKN-16-4-a (Cool wear moisture-wicking fabric)

Lead Free gonad shield SKN-16-5


Lead Free gonad shield SKN-16-5-a

Radiation Protection Accessories

APRON HOLDER



APRON WALL CONSOLE



Apron holder will accommodate up to 2-10 aprons.

Unit is ideal if wall space is limited.

It is designed to prevent declining on side in moving.

Solid stainless steel construction prevents rust in long-term use.

High durability under heavy use.

Can be assembled after delivering.

SK-16-6	up to 10
SK-16-6-a	up to 5
SK-16-6-1	up to 2
SK-16-6-2	up to 5

FOR PREGNANT WOMAN BARRIER



Lead curtain provides separated radiation shielded space by simply drawing.

- protection : 0.5 mmpb
- Order made
- Protection for pregnant woman barrier : 0.5 mmpb
- Order made

SK-16-6	71(L) * 63(W) * 140(H) cm
SK-16-6-a	50(L) * 63(W) * 140(H) cm
SK-16-6-1	20(L) * 50(W) * 25(H) cm
SK-16-6-2	25(L) * 50(W) * 25(H) cm
SK-16-7-1	57(W) * 60(L)/83(H) to 146(H) cm
SK-16-7-3	200(W) * 150(H) cm
SK-16-7-4	100(L) * 65(W) * 70(H) cm
SK-16-7-4-a	100(L) * 65(W) * 70(H) cm

LEAD BARRIER



Lead-glass viewing panel provides exceptional visibility.

- Protection : 1.76 mmpb
- Order made

SK-16-7 lead glass barrier	180(L) * 90(W) cm
Window Size:	80(L) * 90(W) cm
SK-16-7-a lead barrier	180(L) * 90(W) cm
SK-16-7-b	180(L) * 90(W) cm
Window Size:	20(L) * 30(W) cm
SK-16-7-c	180(L) * 90(W) cm
Window Size:	180(L) * 90(W) cm

Radiation Protection Eyewear

Since **human's eyes are more sensitive** to radiation than other parts of human body, **using eyewears or masks is strongly recommended** for users related to the following vocations:
Fluoroscopic Procedures, Cardiac Cath, EP Labs, Urology, Interventional Radiology, Orthopedic Surgery.

- Lead Masks for Full Faces: 0.1 mm Pb
- Lead Eyewears (With Glasses / Without Glasses): 0.75 mm Pb

Lead Masks

SK-400PSM (PANORAMIC FACE SHIELD)



SK-400PSM

The 400PSM Panorama Shield consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from a dual adjustable headpiece.

The acrylic shield can be positioned at any angle to the face so as to minimize scatter radiation from passing under the shield by side adjustable dials. The headgear is adjustable in two positions.

An adjustable headband that saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back. The circumference portion of the headband has an adjustable dial in the rear to increase or decrease the hat size for a comfortable, snug fit.

Wipe clean with a soft cloth. Abrasive materials may scratch front lens.

Cannot be sterilized with high temperatures or high pressure.

- ◆ Weight : 350 grams
- ◆ Protection : .10mm Lead Acrylic
- ◆ Material : Acrylic/Nylon

SK-450FPM (FULL FACE SHIELD)



SK-450FPM

The 450FPM Panorama Mask consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from an adjustable headpiece.

The stationary acrylic shield extends from a headpiece brow connector for a length of 8 inches from top to bottom in front of the entire face. A support strap saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back.

The circumference portion of the headband has an adjustable strap in the rear to increase or decrease the hat size for a comfortable, snug fit. Wipe clean with a soft cloth.

Abrasive materials may scratch front lens. Cannot be sterilized with high temperatures or high pressure.

- ◆ Weight: 450 grams
- ◆ Protection: .10mm Lead Acrylic
- ◆ Material: Acrylic/Nylon

Lead Eyewear (With Glasses)

SK-RE-FGUARD (Beta/Alpha)



The Fitguard fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitguard lead glasses' temples can be benched (adjusted) easily for a secure, comfortable fit. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

The Beta model is ideal for fitting over standard prescription eyewear with square or rounded lenses.

The Alpha model is ideal for fitting over wide prescription eyewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the Fitguard Alpha is our most popular fitover model.

- ◆ Weight: 108g(B) / 91g(A)
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0.50mm Pb Strips
- ◆ Frame Material: Nylon
- ◆ Frame Size: 64-18-150(B) / 69-17-145(A)
- ◆ Color: Black
- ◆ CE: Yes

SK-90 (90 FITOVER)



Designed to fit over prescription eyewear, model 90 fitovers are an excellent choice for those requiring protection against x-ray and liquid splashes.

The large, clear side shields protect the user from lateral exposure while eliminating the "tunnel" effects of look-alikes.

- ◆ Weight: 112g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Standard
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 66-14-120
- ◆ Color: Matte Black
- ◆ CE: Yes

SK-70S (70S ASTRO FLEX) (Without Glasses)



A panoramic goggle constructed of impact resistant propionate.

The frame front has molded side shielding and is fitted to the brow for increased splash protection.

The temple length and frame angle are adjustable for a uniform custom fit.

Model 70S comes with Lead side shields for necessary protection from oblique angle exposure.

- ◆ Weight: 88g
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0.50mm Pb Strips
- ◆ Frame Material: Propionate
- ◆ Frame Size: 54-19-Adj
- ◆ Color: Clear/Blue
- ◆ CE: Yes

Lead Eyewear (Without Glasses)

SK-RE-9941A (9941 ULTRALITE)



The 9941A Ultralites come standard with a rubberized nose bridge. They are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields.

The 9941A is a popular frame and ideal for those who have a high power prescription.

- ◆ Weight: 64g
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0.50mmPb Strips Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 55-22-130
- ◆ Color: Silver
- ◆ CE: Yes

SK-N-BR (NIKE BRAZEN)



The Nike Brazen features an advanced ventilated nose bridge for grip, comfort, and reduced fogging.

Its lightweight nylon design combine streamlined style and enhanced comfort.

The premium Nike Brazen frame offers versatility with its superior coverage and protection.

- ◆ Weight: 75g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 60-16-140
- ◆ Color: Matte Black
- ◆ CE: Yes



SK-RE-RAZER (RAZER)



The Razer's blends a sleek design into a lightweight and durable frame. Its premium features such as a ventilated nose bridge, rubberized temples, and wide-angle lenses maximize comfort and protection from x-ray.

The modern style and use of light-weight materials make the Razer an excellent choice for longprocedures.

Features: Ventilated Nose Bridge, Rubberized Temples for Grip and Comfort, Modern Matte Color Finishes, High Rx Options

- ◆ Weight: 62.5g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb attached
- ◆ Frame Material: Polycarbonate
- ◆ Frame Size: 58-17-135
- ◆ Color: Matte Green, Matte Gray
- ◆ CE: Yes



Non Lead Aprons

Lead Free

Our lighter lead-free material (0.25 / 0.35 / 0.5)

◆ Environmentally Safe

VEST-SKIRT COMBINATION



This Vest-Skirt combination is the preferred choice for complete protection and weight distribution.

The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue.

Skirt features no-slip, quick release belt for added support.

Double padded shoulders for maximum comfort.

Chest pocket standard.

(Pictured with Thyroid Tether.)



◆ Protection : .50mm front / .25mm back

◆ Closures : Velcro and Quick-Release Buckle

VEST-SKIRT COMBINATION

APVS

WRAP AROUND BACK RELIEF



Offering full protection, this Back Relief design transfers the weight from the shoulders to the hips with the aid of a 6" bonded elastic belt.

(See features of our Back Relief belt.)

This apron features Velcro and buckle closures, as well as a Velcro Breakaway side seam for freedom of movement.

Double padded shoulders for added comfort and a chest pocket are standard features.


◆ Protection : .50mm front / .25mm back

◆ Closures : Velcro and Quick-Release Buckle

WRAP AROUND BACK RELIEF

APWABR

BACK RELIEF



Designed to transfer the weight from the shoulders to the hips.

The apron features a unique padded belt system constructed with a bonded elasticized panel.

The stretch proof panel will not lose its shape while providing increased comfort and back support through the life of the apron.

Extra thick shoulder pads, chest pocket and shoulder strap are standard. Optional tether on shoulder strap available upon request.

◆ Protection : .50mm front

◆ Closures : Velcro and Quick Release Buckle

BACK RELIEF APBR

APBR

FLEXBACK



Our Flexback apron is one of the most comfortable models you'll ever wear!

Criss-crossing back panels lift the weight off the shoulders and support the lower back.

Optional shoulder strap available upon request.

Extra thick shoulder pads are standard.

Chest pocket standard.

Shown with Thyroid Collar.

◆ Protection : .50mm front

◆ Closures : Velcro

FLEXBACK

APFB

VEST-SKIRT FLEXBACK



This design combines the benefits of a 2-piece apron with the support of our elasticized flexback apron.

You control the amount of back support by adjusting the Inner elasticized quick release belt.

Chest pocket standard. Shown with Page Cap.

◆ Protection : .50mm front / .25mm back

◆ Closures : Velcro and Quick - Release Buckle

VEST-SKIRT FLEXBACK

APVSFB

BUCKLE APRON



Quick Release Buckle Apron is user friendly with large support straps criss-crossing the back allowing proper placement.

Belt is adjustable in order to accommodate varying sizes. Chest pocket and shoulder pads are standard.

Optional shoulder strap (as shown on Back Relief APBR above) available upon request.

◆ Protection : .50mm front

◆ Closures : Quick Release Buckle

BUCKLE APRON

APBA

FLEXBACK w/SURGICAL DROP SHOULDER



Two great aprons in one! This is the perfect choice for those who want the convenience of a Surgical Drop Apron and the comfort of our Flexback Apron.

Crisscrossing back panels transfer the weight from the shoulders and distribute it to the hips.

Lift the Velcro tabs for quick removal!

Extra thick shoulder pads and chest pocket are standard. Shown with Apron Cap.

◆ Protection : .50mm front

◆ Closures : Velcro

FLEXBACK w/SURGICAL DROP SHOULDER

APFBSB

Non Lead Aprons

SURGICAL DROP-OFF



Fluoro done Remove your apron quickly and easily without breaking the sterile field.

Extra thick shoulder pads for added comfort.

Chest pocket and hanging straps are standard.

◆ Protection : .50mm front

◆ Closures : Velcro

SURGICAL DROP-OFF

APSD

TIE APRON



Our most basic apron. Easy to put on.

Ties are long enough to fit a large spectrum of sizes.

Chest pocket standard.

◆ Protection : .50mm front

◆ Closures : Ties

TIE APRON

APPA

VELCRO ADJUSTABLE



Crisscrossing back straps keep the shoulders from slipping off. Velcro tabs make adjustability quick and easy.

Extra thick shoulder pads, for comfort and chest pocket are standard.

◆ Protection : .50mm front

◆ Closures : Velcro

VELCRO ADJUSTABLE

APVA

MATERNITY APRON



This apron features an adjustable Quick Release Buckle to secure the back panel and easy ties for the front panel allowing complete adjustability during pregnancy. Velcro release shoulders makes removal easy. Apron also includes an additional .50mm equivalent panel.

You can select either the Velcro attachment to the marsupial Inner pouch or Velcro belt.

◆ Protection : .50mm front / 1.0mm(with added panel) and .25 back

◆ Closures : Velcro

MATERNITY APRON

APMA

Lead Gloves and Sleeves

SK-16-1 (Size: 6~9)



Our Soft-Touch finish creates a smooth interior reducing skin irritation. Anatomically correct curved fingers reduce hand fatigue.

Sizing range from 6 to 9 in half sizes for an unequalled comfortable fit. Attenuates up to 50% more direct beam radiation than other leading products. Thinner fingertips provide enhanced tactile sensitivity.

ATTENUATION

RR-2 : 0.012" / 0.30 mm				RR-3 : 0.024" / 0.60 mm			
Kvp	60	90	120	Kvp	60	90	120
Cuff	56%	44%	37%	Cuff	76%	63%	55%
Palm	57%	45%	38%	Palm	75%	62%	54%
Finger	57%	56%	38%	Finger	75%	62%	54%

SKN-16-1 (Size: 6~9)



ATTENUATION

MODEL 1 : 0.009" / 0.22 mm				MODEL 2 : 0.012" / 0.30 mm			
Kvp	60	90	120	Kvp	60	90	120
Cuff	51%	40%	34%	Cuff	62%	51%	43%
Palm	55%	44%	37%	Palm	63%	51%	43%
Finger	55%	44%	37%	Finger	61%	49%	42%

SK-SL (Size: S/M/L)



XGUARD Radiation Reducing Sleeve, when used in conjunction with the XGUARD Radiation Reducing Procedure Glove creates an attenuative barrier from your fingertips to your biceps from scatter radiation.

- ◆ Inner Nitrile coating - reducing proteins, antigens and irritation.
- ◆ Lightweight material reduces arm fatigue.
- ◆ Extreme flexibility does not restrict movement or dexterity.
- ◆ Available in three sizes for a comfortable fit.

ATTENUATION

XGUARD Sleeve	60 Kvp	80 Kvp	100 Kvp	120 Kvp	THICKNESS
	55%	45%	35%	31%	0.012" - 0.30mm


BREAST SHIELD



SURGICAL DROP SHOULDER



THYROID COLLARS




BREAST SHIELD	APSTOLE
SURGICAL DROP SHOULDER	APDRSH
THYROID COLLARS	APTCO / APTCS

Non Lead Aprons Color Swatches

LIMITED EDITION PRINTS*




COOLGUARD PERFORMANCE FABRICS*



SATINGUARD ANTI-MICROBIAL FABRICS*

Featuring stain resistant and anti-microbial technology only available from Protech. Our soft poly material will resist a wide

variety of bacteria for several years.




NYLONS

Unless specified otherwise, nylons come standard 210 denier.



BINDING COLORS



RIPSTOP



Black
#200




Brown
#205




Burgundy
#210




Chocolate
#215




Forest Green
#220



Green
#225



Navy
#230



Purple
#235



Red
#240




Royal Blue
#245




Silver
#250


WIPE-AWAY VINYLs*



Black
#300



Silver
#320




Burgundy Starlight
#305

CLASSIC PRINTS



70's Flowers
#400




Animal Farm
#405




Argyle Green
#410




Asian Garden
#415



Black Lace
#417




Blue Animal River
#420



Blue Flame
#425




Brown Animal Spots
#435



Butterfly
#440



Camo-Back woods
#443



Camo-Cosmo
#450



Camo-Digital Grey
#455



Camo-Pink
#465




Camo-Small Blues
#470



Camo-Small Pinks
#475



Camo-Yellow
#480



Chicklets Blue
#490



Circles
#495



Crayon Hearts
#500



Divebar
#502



Dog Paw
#505



Fire Tongues
#510



Fleur De Lis
#512



Fruity Pebbles
#515




Galaxy-Multi
#520



Galaxy-Royal
#525



Giraffe-Brown
#530




Green Neon Tie-Dye
#540



Hawaiian
#545



Hawaii Postcard
#550



Leopard Print
#560




Lu Lu Grape
#565




Mamba-Blue
#575



Mamba-Pink
#580



Marble Splash
#585




Mardi Gras
#590


CLASSIC PRINTS




Mystery-Multi
#595



Painted Dragon
#600



Paint Splash
#605



Paw Prints
#610



Plaid-Blue
#615



Plaid-Green
#620



Plaid-Light Pink
#625



Plaid-Purple
#630



Plaid-Tan
#635



Plaid-Teal
#640



Purple Neon Tie-Dye
#645




Rainbow
#650



Red Dalmatian
#655



Sheep
#665



Skulls-Gray
#670




Skulls-Red
#675



Snake Skin
#680



Star Gazer
#685



Stars
#690



Steel
#695



Surf Wax
#705




Teal Splash
#710



Tiger
#715



Zebra
#720



Sung Kwang Meditech Co.,Ltd

Information or ordering call

T) 82-31-855-7548 M) 82-10-4590-7548

F) 82-31-855-7549

E) skmeditech@skmeditech.com

Address) 1014-20, Bukhansan-ro, Jangheung-myeon,
Yangju-si, Gyeonggi-do, Korea Sung Kwang Meditech co.,ltd



www.skmeditech.com